

RayDyLy^o

Push Fit cap ready-to-use for vial capping



RayDyLy is used in closed systems, isolators or RABS for **manual** or **automatic** capping as well as for the closing of vials in nests.

Several available options:

- Delivered with or without stopper pre-assembled
- Supplied in bulk or nested
- Gamma or steam sterilized

RayDyLy^o®

The vial capping solution for aseptic filling in isolator



RayDyLy^o Stopper cap

Plastic push fit cap with pre-assembled stopper ready-to-use supplied in bulk or in strips

1



Single use isolator



Discover Disposable Lab

2



Closing 10-caps strip

3



4

RayDyLy^o
more than a product,
a solution



Sterile and depyrogenated vials ready-to-use for aseptic filling

More about RayDyLy^o



ARAYMONDLIFE SASU

1 rue Louis Besançon
38120 Saint-Egrève - FRANCE
Tel +33 4 56 52 61 76
contact@araymond-life.com

www.araymond-life.com
www.raydylyo.com

© *ARAYMONDLIFE SASU - RCS Grenoble 501 397 012 - 1, Rue Louis Besançon - 38120 Saint Egrève - FRANCE has designed this flyer.
** "ARAYMOND Network" means a network of companies which have a license of use of ARAYMOND™. This flyer is provided for information purposes only and does not constitute an offer or an agreement. ARAYMONDLIFE makes no warranty or representation whatsoever, express or implied, including but not limited to the accuracy, reliability, novelty, completeness, fit for a particular purpose or merchantability of the information contained in it. If you need further information, please contact ARAYMONDLIFE.

ARAYMOND 
MORE THAN FASTENING